CHEMOTECHNIQUE
DIAGNOSTICS

BasIQ UltraTM

INSTRUCTIONS FOR USE BasIQ UltraTM Patch Test Units

For professional use only

EN / FR / ES / DE / SE / BG / CZ /
DK GR / EE FI / IE / HR / HU / IT / LT
LV / MT / NO / NL / PL / PT / RO / SK
SI / TR / MY / KR / CN / IL / EG

© 2025 Chemotechnique MB Diagnostics AB
Rev 01/2025

Product description

BasIQ Ultra™ includes 50 Patch Test Units.

BasIQ Ultra™ is intended to be used with topical haptens for the aid of diagnosis of contact allergy.

BasIQ Ultra™ consists of medical grade acrylic adhesive on non-woven carrier tape with 10 chambers made of polyethylene foam coated with medical grade acrylic adhesive with integrated filter paper made of absorbent cellulose paper. The unit is protected by two layers of liners made of siliconized paper.

Indications for use: BasIQ Ultra™ Patch Test Units are intended for use in epicutaneous diagnostic patch testing and clinical research by or under the supervision of a specially trained healthcare professional only. Patch testing is indicated in all cases where contact allergy is suspected or need to be ruled out.

Contraindications: None known.

Warnings and Precautions:

- Apply patch test units on healthy, dry and clean skin.
- Advise patients to refrain from showering, exercise or excessive sweating.
- Advise patients to protect the Patch Test Units from sunlight (UV radiation) during ongoing patch testing.
- Wiping the test site with ethanol before application may improve adhesion of the Patch Test Unit on patients with oily skin. For patients with very dry or oily skin, additional reinforcing tape might be required. Consider clipping (not shaving) the hair on the test site before patch test application.
- Place the Patch Tests on the upper part of the patients back if possible. Avoid applying Patch Test Units to the midline and the scapula. Avoid placing Patch Test Units under a brassiere shoulder band which can cause dislocation of the test units.
- See adverse events section for more warnings.

Special Patient Populations: BasIQ Ultra™ Patch Test Units can be used in all patients consistent with the limitations associated with the use of individual haptens and patch testing in general.

Adverse Events:

- Allergic reactions to tape adhesives are very rare. However, occasional cases of allergic reactions have been reported.
- Irritant reactions caused by the tape adhesive and foam frame adhesive may occur, but usually disappears rapidly.
- Skin may react to the removal of the tests with a slight mechanical irritation, indicated by erythema.
- Itching and burning may occur at the site of application of the Patch Test Units; it can either be caused by

a positive patch test reaction or due to tape irritation.

Instructions for Use (FIG A):

1. Place BasIQ Ultra™ Patch Test Unit on flat surface with the protective liner facing upwards and the arrow pointing forward.
2. Remove the first layer of protective liner to expose the IQ Chambers.
3. Apply hapten in patch test chamber starting with the 1st hapten in the top-right corner. A placement reference guide for is included in the box. The recommended dose is 25 µl*.
4. Gently detach the bottommost part of second layer of liner protecting the backing tape, leaving the top part of the liner still adhered to the tape.
5. Position the patch test by holding it by the top part of the protective liner with the arrow facing upwards and apply the tape onto the patient. Remove the remaining liner completely once the patch test is applied to the patient. Press your palm on the tape for 5 seconds to enhance adhesion.
6. Mark the 1st and the 5th chamber on the skin. Register the test site location in a Patch test record form**.
7. After 48 hours: Remove the Patch Test Units with a gentle diagonal motion to minimize patient discomfort.
8. At removal, in addition to the visual imprints on the skin from each chamber, a slight erythema from the frames should be visible as a sign of good occlusion. The erythema generally subsides within 30 minutes after removal of the test units. The erythema may however prevail for up to 24h or more for patients with sensitive skin. Allow for 15-60 minutes prior to reading initial results, permitting for the resolution of pressure effects.
9. Used Patch Test Units should be folded with the sticky side inwards and be discarded when used in accordance with local regulations.

* 25 µl equals approximately to a 5 mm string of petrolatum or a single drop of liquid.

** Patch Test record forms are available at www.chemotechnique.se

Report all serious incidents related to this product.

Report to the manufacturer at www.chemotechnique.se/products/incident, local distributor and the competent authority in your country.

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